

CASY®cups are measurement cups with lids designed for use with CASY® Cell Counter + Analyzer Systems. The material properties, size and volume of CASY®cups are optimally adapted for cell culture sample preparation. CASY®cups are disposable and designed for single use. Repeated use can have an effect on measurement results.

Description

CASY®cups measure 22 x 65 mm with a filling volume of 15 ml. They are made of Polypropylene (PPN), and are transparent and milky in color. CASY®cups are particle-pure and biologically inert. The material prevents cells from adhering to the cup walls. The cups are heat-resistant up to 80 °C and slowly inflammable.

Packaging

CASY®cups with lids are available directly from innovatis AG in packages of 1,800 cups.

Shelf-life and storage

CASY®cups must be stored at room temperature in a dry, dust-free location. The shelf-life is unlimited.

Waste disposal

Waste must be disposed of in accordance with local cell culture laboratory regulations.

Order number

43003

innovatis AG
CASY®-Technology
Krämerstr. 22
72764 Reutlingen
Germany
Phone: +49 (0)71 21 38786-0
Fax: +49 (0)71 21 38786-99
casy@innovatis.com

innovatis Inc.
301 Lindenwood Drive, Suite 115
Malvern, PA 19355
U.S.A.
Phone: +1-888-283-1564 Toll free
Fax: +1-888-283-1631 Toll free
Phone: +1-610-889-7319
Fax: +1-610-889-7436
usa.info@innovatis.com

www.innovatis.com

Data specification

CASY® cups

cell counter cups with lids



Principle

CASY® supplies are especially designed for use with standardized CASY® equipment. Only the exclusive and regular use of CASY® supplies guarantees long-term consistency of measurement results, as well as long-term stability and optimal functionality of the equipment.



INNOVATIS
QUALITY THAT COUNTS