

CASY®clean is a ready-to-use cleaner for the CASY® Cell Counter + Analyzer Systems.

Description

CASY®clean is a foam-free, alkaline, odorless and colorless cleaning solution for the system and measurement capillaries. CASY®clean is non-toxic and non-flammable.

Properties

Regular use protects the device and measurement capillaries from biological deposits. Through its emulsifying effect, CASY®clean prevents blockages and encrustations of the internal liquid container system and measurement capillaries.

Packaging

CASY®clean is available directly from innovatis AG as a package of 3 x 500 ml bottles.

Shelf-life and storage

The best-before date is indicated on the product packaging. CASY®clean should be stored in a cool, dry place and should be protected from light.

Waste disposal

Waste must be disposed of in accordance with local cell culture laboratory regulations.

Order number

43008

innovatis AG
CASY®-Technology
Krämerstr. 22
72764 Reutlingen
Germany
Phone: +49 (0)71 21 38786-0
Fax: +49 (0)71 21 38786-99
casy@innovatis.com

innovatis Inc.
301 Lindenwood Drive, Suite 115
Malvern, PA 19355
U.S.A.
Phone: +1-888-283-1564 Toll free
Fax: +1-888-283-1631 Toll free
Phone: +1-610-889-7319
Fax: +1-610-889-7436
usa.info@innovatis.com

www.innovatis.com

Data specification

CASY®clean

ready-to-use cleaner for system
and capillaries



Principle

CASY® supplies are especially designed for standardized CASY® equipment. Only the exclusive and regular use of CASY® supplies guarantees long-term consistency of measurement results, as well as long-term stability and optimal functionality of the equipment.



INNOVATIS
QUALITY THAT COUNTS