

CASY[®]ton is a ready-to-use isotonic dilution liquid for cell cultures especially developed for use with CASY[®] Cell Counter + Analyzer Systems. The particle purity of the liquid guarantees a low background to the measurements. CASY[®]ton is manufactured according to the strict standards of the quality management system ISO 9000 ff.

Description

CASY[®]ton is an isotonic, particle-pure, odorless and colorless saline solution.

Properties

CASY[®]ton reduces aggregation of cells in the measuring cup and stabilizes the cells for the duration of the measurement time. The defined conductivity guarantees consistent measurements. Physiological pH and osmolarity ensure the preservation of cell size.

Packaging

CASY[®]ton is available directly from innovatis AG in convenient 10 liter containers, including a plastic dosage tap.

Shelf-life and storage

The best-before date is indicated on the product packaging. CASY[®]ton can be stored at room temperature.

Waste disposal

Waste must be disposed of in accordance with local cell culture laboratory regulations.

Order number

43001

Data specification

CASY[®]ton

ready-to-use dilution liquid
for cell cultures



Principle

CASY[®] supplies are especially designed for standardized CASY[®] equipment. Only the exclusive and regular use of CASY[®] supplies guarantees long-term consistency of measurement results, as well as long-term stability and optimal functionality of the equipment.

innovatis AG
CASY[®]-Technology
Krämerstr. 22
72764 Reutlingen
Germany
Phone: +49 (0)71 21 38786-0
Fax: +49 (0)71 21 38786-99
casy@innovatis.com

innovatis Inc.
301 Lindenwood Drive, Suite 115
Malvern, PA 19355
U.S.A.
Phone: +1-888-283-1564 Toll free
Fax: +1-888-283-1631 Toll free
Phone: +1-610-889-7319
Fax: +1-610-889-7436
usa.info@innovatis.com



INNOVATIS
QUALITY THAT COUNTS

www.innovatis.com