

Data specification

Cedex⁺

Calibration Beads



Product Advantage

- Consistent product quality: produced under GMP conditions
- Beads suspension tested and optimized for carrying out automatic adjustments of the focus and brightness of Cedex and Cedex HiRes systems
- Certificate of Analysis available

Purpose	Beads to calibrate focus and brightness
Size of trading unit	20 mL
Manufacturing	GMP
Product No.	7 60 00 00 16

Product Description

Cedex Calibration beads are used to perform automated focus and brightness adjustments of the Cedex and Cedex HiRes systems. During the adjustment, the microscope stage position is set to ensure the best focus quality and, if necessary, default settings for the illumination intensity are adjusted. The complete routine is described in the Cedex and Cedex HiRes manuals. The procedure takes approximately 15 minutes.

innovatis AG
Meisenstr. 96
33607 Bielefeld
Germany
Phone: +49 (0)521 2997-300
Fax: +49 (0)521 2997-285

innovatis Inc.
301 Lindenwood Drive, Suite 115
Malvern, PA 19355
U.S.A.
Phone: +1-888-283-1564 Toll free
Fax: +1-888-283-1631 Toll free
Phone: +1-610-889-7319
Fax: +1-610-889-7436

info@innovatis.com
www.innovatis.com



INNOVATIS
QUALITY THAT COUNTS