

## Data specification

# Cedex<sup>+</sup>

## Control Beads



### Product Advantage

- Control Beads complement DRSB for the inspection of the suitability of a system.
- Control Beads are cheaper than DRSB and can, therefore, be used considerably more frequently for control measurements (even several times per day).
- The entire linearity range can be checked with Control Beads using three measurement points to confirm that a system is measuring consistently.
- In comparison, one concentration can be checked with DRSB. However, DRSB are very exact and traceable. Multiple systems can, therefore, be adjusted to each other using DRSB (see data specification about DRSB).

<b>Purpose</b>	Beads for frequent three-point density calibration
<b>Density</b>	$5 \times 10^5$ , $1 \times 10^6$ , $5 \times 10^6$ per mL
<b>Size of trading unit</b>	100 mL
<b>Manufacturing</b>	GMP
<b>User manual</b>	yes
<b>Traceability</b>	no
<b>Product No.</b>	7 60 00 00 75 Cedex Control Beads 5E6 $5 \times 10^6$
	7 60 00 00 76 Cedex Control Beads 1E6 $1 \times 10^6$
	7 60 00 00 79 Cedex Control Beads 5E5 $5 \times 10^5$
	7 60 00 00 80 Cedex Control Beads Bundle, incl. all densities: 1 x 100 mL 5E5 1 x 100 mL 1E6 1 x 100 mL 5E6

innovatis AG  
Meisenstr. 96  
33607 Bielefeld  
Germany  
Phone: +49 (0)521 2997-300  
Fax: +49 (0)521 2997-285

7 60 00 00 81  
Cedex Control Beads Bundle,  
incl. all densities:  
3 x 100 mL 5E5  
1 x 100 mL 1E6  
1 x 100 mL 5E6

innovatis Inc.  
301 Lindenwood Drive, Suite 115  
Malvern, PA 19355  
U.S.A.  
Phone: +1-888-283-1564 Toll free  
Fax: +1-888-283-1631 Toll free  
Phone: +1-610-889-7319  
Fax: +1-610-889-7436

info@innovatis.com  
www.innovatis.com

### Product Description

Control Beads are beads that are used for carrying out System Suitability Tests for Cedex and Cedex HiRes systems. They are offered in three concentrations. The exact concentrations are not certified, but are determined for each bottle and Cedex or Cedex HiRes system before use.



**INNOVATIS**  
QUALITY THAT COUNTS