

## Data specification

# Cedex<sup>+</sup>

## Detergent



### Product Advantage

- The kit is ready-to-use and can be installed immediately
- Consistent product quality; produced under GMP conditions
- Time savings: no mixing or stirring required
- Cleaning solution tested and optimized for use in the Cedex and Cedex HiRes systems
- Certificate of Analysis available on request

|                             |                                    |
|-----------------------------|------------------------------------|
| <b>Purpose</b>              | Cleaning of Cedex capillary system |
| <b>Size of trading unit</b> | 1000 mL                            |
| <b>Manufacturing</b>        | GMP                                |
| <b>Product No.</b>          | 7 60 00 00 23                      |

### Product Description

Cedex Detergent is required for cleaning of the capillary system of the Cedex and Cedex HiRes systems. The ready-to-use solution is certified and available as a component of the Reagent Kit and as a separate bottle.

innovatis AG  
Meisenstr. 96  
33607 Bielefeld  
Germany  
Phone: +49 (0)521 2997-300  
Fax: +49 (0)521 2997-285

innovatis Inc.  
301 Lindenwood Drive, Suite 115  
Malvern, PA 19355  
U.S.A.  
Phone: +1-888-283-1564 Toll free  
Fax: +1-888-283-1631 Toll free  
Phone: +1-610-889-7319  
Fax: +1-610-889-7436

info@innovatis.com  
www.innovatis.com



**INNOVATIS**  
QUALITY THAT COUNTS