

## Data specification

# Cedex<sup>+</sup>

## Reagent Kit



### Product Advantage

- The kit is ready-to-use and can be installed immediately (Water must be supplemented!)
- Time savings: No mixing or stirring required
- No effort needed to meet GMP requirements
- Safer handling of hazardous chemicals

#### Cedex Reagent Kit

<b>Purpose</b>	Ready-to-use kit for Cedex Cell Analyzer
<b>Material</b>	Trypan Blue Solution (0,2 %) 2 x 60 mL Detergent (3 % v/v Mucosol) - 60 mL Cleaning Solution - 60 mL Caps for Cedex Reagent Kits Flasks - 3x
<b>Manufacturing</b>	GMP
<b>Product No.</b>	7 60 50 00 00

#### Cedex Refill Packages

<b>Purpose</b>	Refill Packages for Cedex Reagent Kit
<b>Material Product No.</b>	Trypan Blue Solution (0,2 %) – 4 x 60 mL 7 60 50 10 01
<b>Material Product No.</b>	Detergent – 4 x 60 mL 7 60 50 10 02
<b>Material Product No.</b>	Cleaning Solution – 4 x 60 mL 7 60 50 10 03

#### Caps for Cedex Reagent Kit Flasks

If the Reagent Refill Packages and Cedex Reagent Kits are used with the Cedex system, it is necessary to close the bottles with caps containing holes that allow the capillary tubes to pass through in order to avoid contamination.

<b>Purpose</b>	Caps for Cedex Reagent Kit Flasks
<b>Packaging unit</b>	12 Caps
<b>Product No.</b>	7 60 50 10 00

### Product Description

The Cedex Reagent Kit contains all reagents required for preparation, staining and analysis of eukaryotic cells with the Cedex system.

innovatis AG  
Meisenstr. 96  
33607 Bielefeld  
Germany  
Phone: +49 (0)521 2997-300  
Fax: +49 (0)521 2997-285

innovatis Inc.  
301 Lindenwood Drive, Suite 115  
Malvern, PA 19355  
U.S.A.  
Phone: +1-888-283-1564 Toll free  
Fax: +1-888-283-1631 Toll free  
Phone: +1-610-889-7319  
Fax: +1-610-889-7436

info@innovatis.com  
www.innovatis.com



**INNOVATIS**  
QUALITY THAT COUNTS